

# CovaWound™ Foam

## Non-Adherent Foam Dressing

### WHAT IS IT?

CovaWound™ Non-Adherent Foam Dressing is a non-adherent foam dressing with a two layer construction composed of highly absorbent polyurethane foam with a breathable polyurethane film backing. The absorbent polyurethane foam core is designed to absorb and retain exudate, while the breathable outer layer allows moisture vapour transmission away from the dressing.

### HOW DOES IT WORK?

CovaWound™ Non-Adherent Foam Dressing enables healing to progress by creating and maintaining a moist healing environment. The breathable pink outer layer protects against eschar formation and is also waterproof, protecting against external contamination. CovaWound™ Non-Adherent Foam Dressing is soft and conformable providing an optimal moist healing environment and excellent patient comfort.

### WHEN DO I USE IT?

CovaWound™ Non-Adherent Foam Dressing does not adhere to the wound and, therefore, is recommended for wounds surrounded by sensitive skin, including:

- Venous ulcers.
- The management of moderate to heavily exuding wounds such as pressure, leg, diabetic and foot ulcers
- Superficial and partial-thickness burns
- Donor sites
- Trauma wounds
- Post-operative surgical wounds
- CovaWound™ Non-Adherent Foam Dressing can be used under compression bandaging

Please refer to package insert for full instructions.



Reference Number	Size (in/cm)	Dressings per Carton	Cartons per Case	HCPCS
2000505	2in x 2in 5cm x 5cm	10	20	A6209
2000510	2in x 4in 5cm x 10cm	10	20	
2000520	2in x 8in 5cm x 20cm	10	20	
2001010	4in x 4in 10cm x 10cm	10	20	A6209
2001515	6in x 6in 15cm x 15cm	10	10	
2002020	8in x 8in 20cm x 20cm	10	10	

### SUGGESTED WOUND TISSUE TYPE:

NECROTIC (After debridement)	✓
SLOUGHY	✓
GRANULATING	✓
EPITHELIALISING	✓
INFECTED	—

### ABSORPTION LEVEL: HEAVY

