

CovaWound™ Hydrocolloid

Hydrocolloid Dressing



WHAT IS IT?

CovaWound™ Hydrocolloid has a two layer construction composed of an absorbent hydrocolloid matrix laminated to a breathable polyurethane film backing.

HOW DOES IT WORK?

The hydrocolloid matrix absorbs exudate to form a cohesive gel – while the breathable polyurethane film backing provides a high rate of moisture vapour transmission. The combination of these two features ensures that CovaWound™ Hydrocolloid provides an optimal moist wound environment conducive to wound healing and minimizes the chance of damage to healthy peri-wound skin.

WHEN DO I USE IT?

CovaWound™ Hydrocolloid is indicated for the management of light to medium exuding wounds such as:

- Pressure, leg and foot ulcers
- Superficial partial-thickness burns
- Donor sites, diabetic ulcers, postoperative wounds and skin abrasions
- CovaWound™ Hydrocolloid may also be used on small wounds towards the end of the healing phase

Please refer to package insert for full instructions.

Reference Number	Size (in/cm)	Dressings per Carton	Cartons per Case	HCPCS
5000505	2in x 2in 5cm x 5cm	10	40	
5000510	2in x 4in 5cm x 10cm	10	20	
5000520	2in x 8in 5cm x 20cm	10	20	
5001010	4in x 4in 10cm x 10cm	10	20	
5001515	6in x 6in 15cm x 15cm	10	10	
5002020	8in x 8in 20cm x 20cm	10	10	

SUGGESTED WOUND TISSUE TYPE:

NECROTIC	✓
SLOUGHY	✓
GRANULATING	✓
EPITHELIALISING	✓
INFECTED	—

ABSORPTION LEVEL: LOW

