

CovaWound™ SuperAbsorbent

SuperAbsorbent Dressing



WHAT IS IT?

CovaWound™ SuperAbsorbent dressing consists of a soft hydrophilic wound contact layer with a central SuperAbsorbent Polymer (SAP) core which provides excellent fluid absorption and retention.

HOW DOES IT WORK?

The non-adherent hydrophilic wound contact layer enables transmission of fluid through to the SAP core, which locks in fluid even under pressure, decreasing the risk of maceration of the wound and surrounding tissue. The fluid-repellent nonwoven backing acts as a barrier and prevents exudate strike-through.

WHEN DO I USE IT?

CovaWound™ SuperAbsorbent dressing is indicated for the management of moderate to heavily exuding wounds such as:

- Leg ulcers (venous and arterial), foot ulcers, pressure ulcers
- Traumatic and surgical wounds
- First and secondary degree burns
- CovaWound™ SuperAbsorbent dressing can be used in conjunction with compression therapy for venous leg ulcers

Please refer to package insert for full instructions.

Reference Number	Size (in/cm)	Dressings per Carton	Cartons per Case	HPCPS
4007575	3in x 3in 7.5cm x 7.5cm	10	20	
4001010	4in x 4in 10cm x 10cm	10	20	A6251
4001020	4in x 8in 10cm x 20cm	10	20	
4002020	8in x 8in 20cm x 20cm	10	20	A6252
4002030	8in x 12in 20cm x 30cm	10	20	

SUGGESTED WOUND TISSUE TYPE:

NECROTIC (After debridement)	✓
SLOUGHY	✓
GRANULATING	✓
EPITHELIALISING	✓
INFECTED	—

ABSORPTION LEVEL: HEAVY

