

CovaWound™ Foam

Non-Adherent Foam Dressing with Adhesive Border

WHAT IS IT?

CovaWound™ Non-Adherent Foam Dressing with Adhesive Border is a non-adherent foam dressing with a three layer construction composed of highly absorbent polyurethane foam with a breathable polyurethane film, and an adhesive, breathable polyurethane backing. The absorbent polyurethane foam core is designed to absorb and retain exudate, while the breathable outer layers allow moisture vapour transmission away from the dressing.

HOW DOES IT WORK?

CovaWound™ Non-Adherent Foam Dressing with Adhesive Border enables healing to progress by creating and maintaining a moist healing environment. The breathable outer layers protect against eschar formation, and are also waterproof, protecting against external contamination.

CovaWound™ Non-Adherent Foam Dressing with Adhesive Border is soft and conformable providing an optimal moist healing environment and excellent patient comfort.

WHEN DO I USE IT?

CovaWound™ Non-Adherent Foam Dressing with Adhesive Border does not adhere to the wound. It is indicated for the management of light to moderately exuding wounds such as:

- Pressure, leg, diabetic and foot ulcers
- Superficial and partial-thickness burns
- Donor sites
- Trauma wounds
- Post-operative surgical wounds
- CovaWound™ Non-Adherent Foam Dressing with Adhesive Border can be used under compression bandaging

Please refer to package insert for full instructions.



Reference Number	Size (in/cm)	Dressings per Carton	Cartons per Case	HCPCS
2010505	2in x 2in 5cm x 5cm	10	20	
2010510	2in x 4in 5cm x 10cm	10	20	
2010520	2in x 8in 5cm x 20cm	10	20	
2011010	4in x 4in 10cm x 10cm	10	10	

SUGGESTED WOUND TISSUE TYPE:

NECROTIC (After debridement)	✓
SLOUGHY	✓
GRANULATING	✓
EPITHELIALISING	✓
INFECTED	—

ABSORPTION LEVEL: MODERATE

